

# HAMPTON EXTRA ACTIVE MATTRESS



**Product Code:**

HAMHAMP002

**Description:**

- The Hampton Extra is Harvest Healthcare's plus size replacement mattress
- 18 cells made of Nylon with a TPU lining to reduce noise when cells inflate and deflate material
- 4-way stretch PU cover which is vapour permeable and water resistant and has a Crib 5 fire rating
- All cells interchangeable for ease of cleaning
- Foam surround

**Mattress Dimensions:**

2000 x 1200 x 230mm  
2200 x 1200 x 230mm (Also available on request- Comes with Pump 3)

**Weight:**

18kg

**Max. Patient Weight:**

Harvest 3 Pump: 50 Stone (317kg)

**Pressure Relief Risk:**

Very High

**Pump Technical Specification:**

Please see separate pump specification sheet

**Cleaning Guidance:**

Cover should be wiped down with hot water and mild detergent. If heavily soiled, the cover may be machine washed up to 85°C and tumble dried. For disinfection use sodium hypochlorite solution max 1000ppm. Use no bleach, biological or phenolic-based cleaners. Dry completely before storing.

**Warranty**

In terms of Harvest Healthcare Guarantee and Warranties the Harvest Hampton Extra mattress is guaranteed for three years from date of purchase.

**Exclusions:**

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded
3. Components of pumps which are recommended for annual service, are guaranteed for one year.

**Supplier Contact Details:**

t: +44(0)1709 377172      e: sales@harvesthealthcare.co.uk      w: www.harvesthealthcare.co.uk  
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

**Important Safety Notice**

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.  
(Mattress risk categories are to be used as a guide only)

Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

**Version 2: Apr 2018**