

SANDRINGHAM ACTIVE MATTRESS





Product Code:

Description:

- The Sandringham is Harvest Healthcare's replacement system with cells and foam
- 18 cells made of TPU material reducing noise while the cells inflate and deflate, on 50mm of foam
- Pillow section for patient stability
- All cells are interchangeable for ease of repair and cleaning
- Light mattress allowing for the carer to easily move the mattress
- 4-way stretch welded PU cover which is vapour permeable and water resistant and has a Crib 5 fire rating
- Compatible with all of Harvest Healthcare's pumps allowing for a range of maximum patient weights

Mattress Dimensions: 2000 x 900 x 160mm

Weight: 9kg

Max. Patient Weight: Harvest 1 Pump: 25 Stone (158 kg)

Harvest 250 Pump: 31 Stone (200kg)

Pressure Relief Risk: High

Pump Technical Specification: Please see separate pump specification sheet

Cleaning Guidance:

Cover should be wiped down with hot water and mild detergent. If heavily soiled, the cover may be machine washed up to 65°C and tumble dried. For disinfection use sodium hypochlorite solution max 1000ppm. Use no bleach, biological or phenolic-based cleaners. Dry completely before storing.

Warranty

In terms of Harvest Healthcare Guarantee and Warranties the Harvest Sandringham mattress is guaranteed for three years from date of purchase.

Exclusions:

- 1. Removable covers are guaranteed for one year
- 2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded

Supplier Contact Details:

t: +44(0)1709 377172 e: sales@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60